

## List of Old Models

Model numbers beginning with “S” or “XL” were manufactured under the Misonix brand. All other models were made under the Heat Systems Ultrasonics brand. Over the years we have private labeled sonicators for Fisher Scientific and Virtis. Contact us if you have any questions.

<u>Model #</u>	<u>Power/Watts</u>	<u>Features/Converter</u>	<u>Years Made</u>
S-3000	600	Auto-Tuning/CL5	2000-Present
XL-2020	550	Pulsing & Timing/CL4	1990-2000
W-385	374	Pulsing & Timing/CL3	1986-1990
W-375	350	Pulsing & Timing/CL2	1983-1985
W-225	200	Pulsing & Timing/H-1, CL2	1978-1982
XL-2015	475	Pulsar Only / CL4	1990-1995
XL-2010	475	Manual Oper / CL4	1990-1993
W-380	375	Pulsar Only / CL4	1986-1990
W-370	350	Timer Only / CL2	1978-1985
W-220	200	Timer Only / H-1, CL2	1978-1985
XL-2000	100	Pulse/Continuous/CML-4	1997-Present
XL-2007	100	Manual F/S/CML-3	1995-1997
XL2005	50	Manual Pulse / CML-2	1990-1995
MS-50	50	Manual Pulse / CML-2	1985-1990
MS-25	25	Manual Pulse / CML-2	1984-1985
MS-10	10	Manual Pulse / M-R	1984-1985
LS-75	75	Manual Pulse / M-R	1971-1984
W-350	350	Pulsing Only / H-1A	1972-1978
W-200 P	150	Pulsing Only / “L”	1970-1976
W-200 R	200	Pulsing Only / H-1	1976-1978
W-185 C,D,E	150	Manual Pulse / “L”	1968-1970
W-185 F	150	Manual Pulse / H-1	1976-1978
W-140 C,D,E	100	Manual Pulse / “L”	1965-1968

Units with an “F” prefix such as F-550 and F-60 were made by Misonix but sold under the Fisher Scientific brand name.

Qsonica can provide replacement probes for most units. Repair services are available for most items made after 1990. Trade-in discounts are available if you wish to upgrade your old unit for a newer model.